


## 自动回缩安全注射器 Auto retractable safety syringe

### + 产品优势 / Product advantages

- 单手安全操作，单手激活
- 药物注射完成后自动回缩
- 自动回缩后保护针尖不暴露
- 无需增加额外操作步骤，减少培训成本
- 固定针头，无特别死腔，药液残留量低
- 减少废弃物体积及处理成本
- Single use safety with one handed operation;
- Fully auto retraction after medication is discharged;
- Non-exposure of needle after automatic retraction;
- Requires minimum training;
- Fixed needle, no dead space;
- Reduces disposal size and cost of waste disposal;



Also available with other needle size on request  
针头也可根据要求更换大小

规格 / Specification				
0.5ml	1ml	2-3ml	5ml	10ml

针头规格 / Needle gauge					
18G	19G	20G	21G	22G	23G
24G	25G	26G	27G	29G	30G

### + 操作参考图 / How to use

<p>按标准程序抽取药液 Draw up medicine using standard procedure</p> 	<p>使用无菌技术进行注射 Give injection using aseptic techniques</p> 
<p>完全按压芯杆激活回缩机制 Fully depress the plunger to activate retraction</p> <p>针自动缩回到针筒并锁定 Needle automatically retract into plunger &amp; locked</p> 	<p>将注射器丢弃进锐器箱 Dispose of syringe in sharps container</p> 